
lnBody570

The Fastest, Upgraded Solution for Monitoring Your Health



See What You're Made of

Monitoring weight is not enough to see progressive changes in health and body



* Height: 174cm, Age: 27, Gender: Male

Weight alone does not correctly reflect the effects of exercise and improved diet. The graph above shows a man whose weight, throughout a month of exercise, had minimal change; however, his muscle mass and fat mass dramatically increased and decreased respectively.

Changes in muscle and fat mass are vital to understanding the body's true composition. Upon this, the InBody Test utilizes its patented technologies to reveal segmental body fat and muscle distribution percentages as well as the body water balance.

The InBody Test shows a true assessment of the body.

InBody, the Body Composition Analyzer

Have an effective exercise plan and track the progress of the body's change with the Inbody Test

- · More than 40 result outputs are given through an easy and fast InBody Test.
- The InBody results are used as the first screening tool for indicators of potential diseases and poor health.
- · Segmental Muscle Analysis allows for a more focused exercise plan.
- · Body Water Analysis can be an indicator of a poor physical status.



Accuracy and Reliability of the InBody are Proven by the World's Top Journals and Scholars

More than 500 articles have been published by renowned journals

Clinical reliability was proved by the world's medical professionals in numerous articles.

The InBody has 98% of correlation with the gold standard device DXA and the InBody's own technologies hold patents in numerous countries throughout the world.

The InBody Technology

Arms, trunk, and legs are measured separately

High precision by using a set of high and low frequencies simultaneously

Highly reproducible data due to fixed measuring locations on the wrist and ankle



No need of empirical estimation

Age or gender does not affect the result

The InBody's body composition data deliver research-level results and thus have been utilized by thousands of studies to accurately track changes in body composition.

Validation Studies

Kriemler, S., Puder, J., Zahner, L., Roth, R., Braun-Fahrländer, C., & Bedogni, G. (2008). Cross-validation of bioelectrical impedance analysis for the assessment of body composition in a representative sample of 6-to 13-year-old children. *European journal of clinical nutrition*, 63(5), 619-626.

Ling, C. H., de Craen, A. J., Slagboom, P. E., Gunn, D. A., Stokkel, M. P., Westendorp, R. G., & Maier, A. B. (2011). Accuracy of direct segmental multi-frequency bioimpedance analysis in the assessment of total body and segmental body composition in middle-aged adult population. *Clinical Nutrition*, 30(5), 610-615.

Lim, J. S., Hwang, J. S., Lee, J. A., Kim, D. H., Park, K. D., Jeong, J. S., & Cheon, G. J. (2009). Cross-calibration of multi-frequency bioelectrical impedance analysis with eight-point tactile electrodes and dual-energy X-ray absorptiometry for assessment of body composition in healthy children aged 6–18 years. *Pediatrics International*, 51(2), 263-268.

Utter, A. C., & Lambeth, P. G. (2010). Evaluation of multifrequency bioelectrical impedance analysis in assessing body composition of wrestlers. *Med Sci Sports Exerc*, 42(2), 361-7.


The InBody570, Your One and Only Smart Healthcare Solution

Just step on and let the InBody570 do the rest




Fast and easy test

User friendly interface with voice guidance lets you easily take the InBody Test and collect results.



Two different test modes: Self Mode and Professional Mode

Two different modes satisfy both the user and the consultant. The user can easily take the test with the Self Mode, by only inputting their own height. When the Professional Mode is on, a more detailed consultation information is provided on the screen.




The InBody Results Sheet with more than 40 outputs

Various body composition outputs are provided on the single-paged InBody Results Sheet. Provide individualized consultation by customizing outputs on the InBody Results Sheet and track progress with the Body Composition History graph. The personalized InBody Results Sheet will give enough motivation to exercise!

Lookin'Body Data Management Software The Best Way to Manage from Your PC


Wireless connection between the InBody570 and PC allows for better data management



Wireless Connection with the InBody570

Connect your PC with the InBody570 via Wi-Fi or Bluetooth. User data will be listed up in your PC and by using it, you can remotely control the InBody570, save details of the user, and manage appointments with email service.

* Lookin'Body is an optional software.



Strategic Consultation

The Body Composition History graph of each category helps you see your body composition change at a glance. Additionally, the comment functionality of each consultation allows for a more personalized healthcare.



InBody

[InBody570]

BIOSPACE

TEL:02-501-3939 FAX:02-501-2716

1 Body Composition Analysis

	Values	Total Body Water	Soft Lean Mass	Fat Free Mass	Weight
Total Body Water (L)	27.5 $(26.3 \sim 32.1)$	27.5	35.1		
Protein (kg)	7.2 (7.0 ~ 8.6)		$(33.3 \sim 40.7)$	37.3 (35.8 ~ 43.7)	59.1
Minerals (kg)	2.63 (2.44 ~ 2.98)	non-osseous		(33.6 ~ 43.7)	(43.9 ~ 59.5)
Body Fat Mass (kg)	21.8 (10.3 ~ 16.5)				

2 Muscle-Fat Analysis

		Uı	nder		Vorma				Ov	er			
Weight	(kg)	55	70	85	100	115 5 9	.130	145	160	175	190	205	96
SMM Skeletal Muscle Mass	(kg)	70	80	90 19	9.6	110	120	130	140	150	160	170	96
Body Fat Mass	(kg)	40	60	80	100	160	■ 21.8	280	340	400	460	520	96

3 Obesity Analysis

		U	nder		Norma	ıl 📗			Ov	er		
BMI Body Mass Index	(kg/m²)	10.0	15.0	18.5	21.0	^{25.0} 24	30.0	35.0	40.0	45.0	50.0	55.0
PBF Percent Body Fat	(%)	8.0	13.0	18.0	23.0	28.0	33.0	38.0 ■ 36.9	43.0	48.0	53.0	58.0

4 Segmental Lean Analysis Based on ideal weight Based on current weight

8		U	nder		Norma	1		Ü	Ov	er		Ŭ	
Right Arm	(kg) (%)	40	60	80	100	2.02	140	160	180	200	220	240	96
Left Arm	(kg) (%)	40	60	80	100 1. 98.		140	160	180	200	220	240	96
Trunk	(kg) (%)	70	80	90	100 17. 95.4	7 110	120	130	140	150	160	170	96
Right Leg	(kg) (%)	70	80	= 5.20 83.6	0 100	110	120	130	140	150	160	170	96
Left Leg	(kg) (%)	70		5.02 0.6	100	110	120	130	140	150	160	170	96

5 Body Water Analysis

	U	nder		Vorma	ı			Ov	er		
	0.320	0.340	0.360	0.380	0.390	0.400	0.410	0.420	0.430	0.440	0.450
ECW Ratio						0.39	97				

6 Body Composition History

Weight (kg)	65.3	63.9	62.4	61.8	62.3	60.9	60.5	59.1
SMM Skeletal Muscle Mass (kg)	20.1	20.0	19.7	19.7	19.8	19.7	19.8	19.6
PBF Percent Body Fat (%)	41.3	40.7	39.2	39.0	39.4	38.6	37.8	36.9
ECW Ratio	0.399	0.398	0.396	0.396	0.397	0.396	0.398	0.397
▼ Recent □ Total	11.10.10 09:15	11.10.30 09:40	11.11.02 09:35	11.12.15 11:01	12.01.12 08:33	12.02.10 15:50	12.03.15 08:35	12.05.04 09:46

InBody Score—

68/100 Points

* Total score that reflects the evaluation of body composition. A muscular person may score over 100 points.

Weight Control

Target Weight	51.7 kg
Weight Control	- 7.4 kg
Fat Control	- 9.9 kg
Muscle Control	+2.5 kg

Obesity Evaluation—

В	M	I	□ Under	Mormal	□ Over □ Extremely Over
					Over

P B F □ Normal □ Over ★ Extremely Over

Body Balance Evaluation —

Upper	Balanced Slightly Extremely Unbalanced
Lower	☐ Balanced Slightly ☐ Extremely Unbalanced
Upper- Lower	☐ Balanced Slightly ☐ Extremely Unbalanced

Segmental Fat Analysis —

	V — A
Right Arm	(1.5kg)————————————————————————————————————
Left Arm	(1.6kg)————————————————————————————————————
Trunk	(11.7kg)————240.0%
Right Leg	(2.9kg)———132.0%
Left Leg	(2.9kg) = 132.0%

Additional Data

-			
	Intracellular Water	16.6 L	(16.3 ~ 19.9)
	Extracellular Water	10.9 L	$(10.0 \sim 12.2)$
	Basal Metabolic Rate	1176 kcal	
	Waist-Hip Ratio	0.92	$(0.75 \sim 0.85)$
	Visceral Fat Level	12	(1~9)
	Obesity Degree	114 %	(90~110)
	Bone Mineral Content	2.18 kg	$(2.01 \sim 2.45)$
	Body Cell Mass	23.8 kg	$(23.4 \sim 28.6)$
	Arm Circumference	30.2 cm	
	Arm Muscle Circumference	25.7 cm	

Results Interpretation QR Code

Scan the QR Code to see results interpretation in more detail.



(B) Impedance —

		RA	LA	TR	RL	LL
$\mathbf{Z}(\Omega)$	5kHz	373.1	385.4	25.7	303.0	314.1
	50 kHz 50 kHz 500 kHz	337.2	352.5	23.0	282.3	289.8
	500kHz	297.4	311.5	19.1	258.1	267.8

The InBody Results Sheet

Body composition assessment and nutritional information at a glance

1 Body Composition Analysis

The body weight is the sum of Total Body Water, Protein, Minerals and Body Fat Mass. Maintain a balanced body composition to stay healthy.

2 Muscle-Fat Analysis

Compare the bar lengths of Skeletal Muscle Mass and Body Fat Mass. The longer the Skeletal Muscle Mass bar is compared to the Body Fat Mass bar, the stronger the body is.

3 Obesity Analysis

BMI is an index used to determine obesity by using height and weight. PBF is the percentage of body fat compared to body weight.

4 Segmental Lean Analysis

Evaluates whether the muscles are adequately developed in the body. The top bar shows the comparison of the muscle mass to the ideal weight while the bottom bar shows that of the current weight.

5 Body Water Analysis

ECW Ratio, the ratio of Extracellular Water to Total Body Water, is an important indicator whether the body water is balanced.

6 Body Composition History

Track the history of the body compositional change. Take the InBody Test periodically and monitor the progress.

| Indices | Page | Page

1 InBody Score

This score shows the evaluation of the body composition, which includes muscle, fat, and water in the body.

8 Weight Control

See how the body measures up to the recommended Weight, Muscle Mass, and Body Fat Mass for a good balance. The '+' means to gain and the '-' means to lose.

9 Obesity Evaluation

Evaluates obesity based on your BMI and Percent Body Fat.

10 Body Balance Evaluation

Evaluates the balance of the body based on the Segmental Lean Analysis.

11 Segmental Fat Analysis

Evaluates whether the amount of fat is adequately distributed in all parts of the body. Each bar shows fat mass in comparison to the ideal.

Additional Data

Various nutritional outputs are provided such as Intracellular Water, Extracellular Water, Basal Metabolic Rate, Waist-Hip Ratio, Visceral Fat Level, Obesity Degree, and so on.

13 Impedance

Impedance is the resistance value measured when electrical currents are applied throughout the body. Based on the measured data, key body composition outputs can be analyzed. Impedance is also used for many research purposes.

The InBody Results Sheet for a Child

Specially designed results sheet with Growth Graph is available for a Child

InBody570 Specifications

Key Specifications

Bioelectrical Impedance Analysis (BIA) Measurement Items

15 Impedance Measurements by Using 3 Different Frequencies (5kHz, 50kHz, 500kHz) Bioelectrical Impedance (Z) at Each of 5 Segments (Right Arm, Left Arm, Trunk, Right Leg, and Left Leg)

Electrode Method Tetrapolar 8-Point Tactile Electrodes

Measurement Method Direct Segmental Measurement Bioelectrical Impedance Analysis Method (DSM-BIA)

Body Composition Calculation Method Outputs

(InBody Results Sheet)

No Empirical Estimation Results and Interpretations

Body Composition Analysis (Total Body Water, Protein, Soft Lean Mass, Minerals, Body Fat Mass, Weight), Muscle-Fat Analysis (Weight, Skeletal Muscle Mass, Body Fat Mass), Obesity Analysis (Body Mass Index, Percent Body Fat), Segmental Lean Analysis (Based on ideal weight/Based on current weight; Right Arm, Left Arm, Trunk, Right Leg, Left Leg), Body Water Analysis (ECW Ratio), Body Composition History (Weight, Skeletal Muscle Mass, Percent Body Fat, ECW Ratio), InBody Score, Body Type (Based on BMI and Percent Body Fat, Athletic Shape, Slightly Obese, Obesity, Muscular Shape, Average, Slightly Obese, Slim Muscular, Slim Sarcopenic Obesity, Thin, Slightly thin), Weight Control (Target Weight, Weight Control, Fat Control, Muscle Control), Obesity Evaluation (BMI, Percent Body Fat), Nutrition Evaluation (Protein, Minerals, Fat Mass), Body Balance (Upper, Lower, Upper-Lower), Segmental Fat Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), Segmental Extrados (Neck, Chest, Abdomen, Right Arm, Left Arm, Hip, Right Thigh, Left Thigh), Additional Data (Intracellular Water, Extracellular Water, Skeletal Muscle Mass, Basal Metabolic Rate, Waist-Hip Ratio, Waist Circumference, Visceral Fat Level, Obesity Degree, Bone Mineral Content, Body Cell Mass, Arm Circumference, Arm Muscle Circumference), Blood Pressure (Systolic, Diastolic, Pulse,

Mean Artery Pressure, Pulse Pressure, Rate Pressure Product)

· Results Interpretation QR Code

·Impedance

(InBody Results Sheet for a Child)

Results and Interpretations

Body Composition Analysis (Total Body Water, Protein, Soft Lean Mass, Minerals, Body Fat Mass, Weight), Muscle-Fat Analysis (Weight, Skeletal Muscle Mass, Body Fat Mass), Obesity Analysis (Body Mass Index, Percent Body Fat), Growth Graph (Height, Weight), Body Composition History (Height, Weight, Skeletal Muscle Mass, Percent Body Fat), Growth Score, Obesity Evaluation (BMI, Percent Body Fat), Nutrition Evaluation (Protein, Minerals, Fat Mass), Body Balance (Upper, Lower, Upper-Lower), Segmental Lean Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), Additional Data (Intracellular Water, Extracellular Water, Basal Metabolic Rate, Child Obesity Degree, Bone Mineral Content, Body Cell Mass), Blood Pressure (Systolic, Diastolic, Pulse, Mean Artery Pressure, Pulse Pressure, Rate Pressure Product)

· Results Interpretation QR Code

·Impedance

Feature Specifications

Optional Equipment Stadiometer from BIOSPACE and blood pressure monitors from BIOSPACE Name, Address, and Contact Information can be shown on the InBody Results Sheet Custom Logo

Digital Results LCD Monitor, Data management software Lookin'Body Types of Result Sheets InBody Test Results Sheet, InBody Test Results Sheet for a Child

Voice Guidance Provides audible indication for test in progress, test complete, and successfully saved settings changes Database Test Results can be saved when inputting ID into the InBody. Test Results can be saved up to 100,000

Self Mode, Professional Mode Test Mode

Settings: Configure settings and manage data Administrator Menu

Troubleshooting: Additional information to help use the InBody570

USB Thumb Drive Copy, backup, or restore the InBody570 data (data can be viewed on Excel or Lookin'Body data management software) Backup data saved in the InBody by using a USB Thumb Drive, Restore results on the InBody from a backup file Backup Data

Other Specifications

Applied Rating Current $400 \mu A (\pm 40 \mu A)$

Adapter Manufacture BridgePower Corp. JMW140KA1240F02 Model

AC 100 ~ 240V, 50/60Hz, 1.2A Power Input

Power Output DC 12V, 3.4A 800 × 480 7inch Color TFT LCD

Display Type Internal Interface Touchscreen, Keypad

External Interface RS-232C 4EA, USB HOST 2EA, USB SLAVE 1EA, LAN (10T) 1EA, Bluetooth 1EA, Wi-Fi 1EA

Compatible Printer Laser/Inkjet Printers (Printers recommended by BIOSPACE)

* A list of printers compatible with the InBody 570 can be found at http://www.inbodyservice.com

Dimension 522 (W) × 893 (L) × 1113 (H): mm

20.55 (W) × 35.16 (L) × 43.82 (H): inch

Equipment Weight 24kg (52.9lbs) Testing Time About 50 seconds

Operation Environment $10 \sim 40$ °C (50 ~ 104°F), $30 \sim 75$ % RH, $70 \sim 106$ kPa

-20 ~ 70°C (-4 ~ 158°F), 10 ~ 95% RH, 50 ~ 106kPa (No Condensation) Storage Environment

Testing Weight Range $10 \sim 250 \text{kg} \ (22.0 \sim 551.11 \text{bs})$

Testing Age Range 3 ~ 99 years

Height Range 95 ~ 220cm (3ft. 1.40in. ~ 7ft. 2.61in.)

* Specifications may change without prior notice.

BIOSPACE is a body composition analysis device manufacturer that has acquired over 80 patent rights across the globe.









522







893



BIOSPACE

Biospace Co., Ltd. [HEAD OFFICE] Biospace, Inc. [USA] TEL: +82-2-501-3939 TEL: +1-323-932-6503 FAX: +82-2-578-2716 FAX: +1-323-952-5009

Website: http://www.e-inbody.com Website: http://www.biospaceamerica.com E-mail: info@inbody.com E-mail: USA@biospaceamerica.com

Biospace Japan Inc. [JAPAN]

TEL: +81-03-5298-7667 FAX: +81-03-5298-7668 Website: http://www.inbody.co.jp E-mail: inbody@inbody.co.jp

Biospace China. [CHINA]

TEL: +86-21-64439738, 9739, 9705

FAX: +86-21-64439706

Website: http://www.biospacechina.com E-mail: info@biospacechina.com